**Baxter**

# HDx THERAPIE

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<sup>++</sup> Angepasst nach Blackowicz: Economic evaluation of expanded hemodialysis with the **THERANOVA 400** dialyzer: A post hoc evaluation of a randomized clinical trial in the United States.<sup>31</sup>


<sup>+++</sup> Angepasst nach Ariza: An initial evaluation of **HDx** on hospitalizations, drug utilization, costs, and patient utility in Colombia.<sup>16</sup>

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Angepasst nach Rosner M, et al. Classification of Uremic Toxins and Their Role in Kidney Failure. Clin J Am Soc Nephrol. 2021;16(12):1918-1928<sup>32</sup>  
EUTOX Uremic Solutes Database. June 2022. Uremic-toxins.org<sup>50</sup>

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## REFERENZEN

1. Zweigart C, Boschetti-de-Fierro A, Hulko H, et al. Medium cut-off membranes – closer to the natural kidney removal function. *Int J Artif Organs*. 2017; 40(7):328–334.
2. Ronco, C. The Rise of Expanded Hemodialysis. *Blood Purif*. 2017;44:1-VIII.
3. Molano-Trivino A, Sanabria M, Vesga J, Buitrago G, Sanchez R, Rivera A. Effectiveness of medium cut-off vs high flux dialyzers: a propensity score matching cohort study. *In Nephrol Dial Transport*. 2021;36:486-U948.
4. Penny JD, Jarosz P, Salerno FR, Lemoine S, McIntyre CW. Impact of Expanded Hemodialysis Using Medium Cut-off Dialyzer on Quality of Life: Application of Dynamic Patient-Reported Outcome Measurement Tool. *Kidney Medicine*. 2021;3(6):992-1002.
5. Hutchison CA, Wolley M. The Rationale for Expanded Hemodialysis Therapy (HDx). *Contrib Nephrol*. 2017; 191:142-52.
6. Kirsch AH, Lyko R, Nilsson LG, et al. Performance of hemodialysis with novel medium cut-off dialyzers. *Nephrol Dial Transpl*. 2017; 32(1):165-72.
7. Zhang JC, El-Majzoub S, Li M, et al. Could symptom burden predict subsequent healthcare use in patients with end stage kidney disease on hemodialysis care? A prospective, preliminary study. *Ren Fail*. 2020;42(1):294-301.
8. Lim JH, Park Y, Yook JM, et al. Randomized controlled trial of medium cut-off versus high-flux dialyzers on quality of life outcomes in maintenance hemodialysis patients. *Sci Rep*. 2020;10(1):1-11.
9. Alarcon J.C, Bunch A, Ardila F, et al. Impact of Medium Cut-Off Dialyzers on Patient-Reported Outcomes: COREXH Registry. *Blood Purif*. 2021;50:110-118.
10. Anwar N, Naz N, Reynolds A, Chamber. HDx: is it a better way to dialyse? *Nephrol Dial Transplant*. 2020;35(3):1109.
11. Bolton S, Gair R, Nilsson LG, Matthews M, Stewart L, McCullagh N. Clinical Assessment of Dialysis Recovery Time and Symptom Burden: Impact of Switching Hemodialysis Therapy Mode. *Patient Relat Outcome Meas*. 2021;12:315-321.
12. Rayner HC, Zepel L, Fuller DS, et al. Recovery time, quality of life, and mortality in hemodialysis patients: the Dialysis Outcomes and Practice Patterns Study (DOPPS). *Am J Kidney Dis*. 2014;64:86–94.
13. Cozzolino M, Magagnoli L, Ciceri P, Conte F, Galassi A. Effects of a medium cut-off (ThERANOVA®) dialyser on haemodialysis patients: a prospective, cross-over study. *Clin. Kidney J*. 2021;14(1):382-389.
14. Sanabria RM, Hutchison CA, Vesga, JI, Ariza JG, Sanchez R, Suarez AM. Expanded Hemodialysis and Its Effects on Hospitalizations and Medication Usage: A Cohort Study. *Nephron*. 2021;145(2):179-187.
15. Lim JH, Jeon Y, Yook JM, et al. Medium cut-off dialyzer improves erythropoiesis stimulating agent resistance in a hepcidin-independent manner in maintenance hemodialysis patients: results from a randomized controlled trial. *Sci Rep*. 2020;10(1):1-10.
16. Ariza JG, Walton SM, Suarez AM, Sanabria M, Vesga JI. An initial evaluation of expanded hemodialysis on hospitalizations, drug utilization, costs, and patient utility in Colombia. *Ther Apher Dial*. 2021;25(5):621-627.
17. Hadad-Arrascue F, Nilsson LG, Rivera AS, Bernardo AA, Cabezuero Romero JB. Expanded hemodialysis as effective alternative to on-line hemodiafiltration: A randomized mid-term clinical trial. *Ther Apher Dial*. 2022;26(1):37-44.
18. Pisoni RL, Wikström B, Elder SJ, et al. Pruritus in haemodialysis patients: International results from the Dialysis Outcomes and Practice Patterns Study (DOPPS). *Nephrol Dial Transplant*. 2006;21(12):3495-3505.
19. Su H, Lei C.T, Zhang C. Interleukin-6 Signaling Pathway and Its Role in Kidney Disease: An Update. *Front Immunol*. 2017;8:405.
20. Kimmel, M, Alscher DM, Dunst R, Braun, et al. The role of micro-inflammation in the pathogenesis of uraemic pruritus in haemodialysis patients. *Nephrology Dialysis Transplantation*. 2006;21(3):749-755.
21. Lin XW, Zhang JF, Qiu MY, et al. Restless legs syndrome in end stage renal disease patients undergoing hemodialysis. *BMC Neurol*. 2019:19:47.
22. Freire de Menezes A, et al. *Int J Nephrol*. 2018;25:2018:1414568.
23. Giannaki CD, et al. Epidemiology, impact and treatment options of restless legs syndrome in end-stage renal disease patients; an evidence based review. *Kidney Int*. 2014;85(6):1275-1282.
24. La Manna G, Pizza F, Persici E et al. Restless legs syndrome enhances cardiovascular risk and mortality in patients with end-stage kidney disease undergoing long- term haemodialysis treatment. *Nephrol Dial Transplant*. 2011;26:1976–1983.
25. Olsson M.G, Allhorn M, Bulow L, et al. Pathological conditions involving extracellular hemoglobin molecular mechanisms, clinical significance, and novel therapeutic opportunities for alpha-1-microglobulin. *Antioxid Redox Signal*. 2012;17:813–846.
26. Bossola M, Tazza L. Postdialysis fatigue: a frequent and debilitating symptom. *Semin Dial*. 2016;29:222–227.
27. Bossola M, Vulpio C, Tazza L. Fatigue in chronic dialysis patients. *Semin Dial*. 2011;24:550–555.
28. Lindsay RM, Heidenheim PA, Nesrallah G, Garg AX, Suri R. Daily hemodialysis study group London health sciences centre: Minutes to recovery after a hemodialysis session: a simple health-related quality of life question that is reliable, valid, and sensitive to change. *Clin J Am Soc Nephrol*. 2006;1:952–959.
29. Bossola M, Di Stasio E, Monteburini T, et al. Recovery Time after Hemodialysis Is Inversely Associated with the Ultrafiltration Rate. *Blood Purif*. 2019;47(1-3):45-51.
30. Bossola M, Di Stasio E, Giungi S, Rosa F, Tazza L. Fatigue is associated with serum interleukin-6 levels and symptoms of depression in patients on chronic hemodialysis. *J Pain Symptom Manag*. 2015;49:578–585.
31. Blackowicz MJ, Falzon L, Beck W, Tran H, Weiner DE. Economic evaluation of expanded hemodialysis with the ThERANOVA 400 dialyzer: A post hoc evaluation of a randomized clinical trial in the United States. *Hemodialysis International*. 2022. <https://doi.org/10.1111/hdi.13015>.
32. Rosner M, Reis T, Husain-Syed, et al. Classification of Uremic Toxins and Their Role in Kidney Failure. *Clin J Am Soc Nephrol*. 2021;16(12):1918-1928.
33. Ronco C, Marchionna N, Brendolan A, Neri M, Lorenzin A, Martínez Rueda AJ. Expanded haemodialysis: from operational mechanism to clinical results. *Nephrology Dialysis Transplantation*. 2018;33(3):iii41–iii47.
34. Wolley M, Jardine M, Hutchison CA. Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol*. 2018;13(5):805-814.
35. Ronco,C. *Expanded Hemodialysis: Innovative Clinical Approach in Dialysis*. Vol.191. Karger Medical and Scientific Publishers. 2017.
36. Sakurai K. Biomarkers for Evaluation of Clinical Outcomes of Hemodiafiltration. *Blood Purif*. 2013;35(1):64-68.
37. Lorenz G, Schmalenberg M, Kemmner S, et al. Mortality prediction in stable hemodialysis patients is refined by YKL-40, a 40-kDa glycoprotein associated with inflammation. *Kidney Int*. 2018;93(1):221-230.
38. Boschetti-de-Fierro A, Voigt M, Storr M, et al. MCO Membranes: Enhanced Selectivity in High-Flux Class. *Sci Rep*. 2015;5(1):18448.
39. Ronco C, Clark WR. Haemodialysis membranes. *Nat Rev Nephrol*. 2018;14(6):394-410.
40. Device Classification Under Section 513(f)(2)(De Novo). <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/denovo.cfm?id=DEN190042>. Published 2020. Accessed April 26, 2022.
41. Mazairac AH, Blankestijn PJ, Grooteman MP, et al. The cost-utility of haemodiafiltration versus haemodialysis in the Convective Transport Study. *Nephrol Dial Transplant*. 2013;28(7):1865-1873.
42. Schepers E, Glorieux G, Eloot S, et al. Assessment of the Association Between Increasing Membrane Pore Size and Endotoxin Permeability Using a Novel Experimental Dialysis Simulation Set-Up. *BMC Nephrol*. 2018;19:1.
43. Boschetti-de-Fierro A, Beck W, Krause B, Hildwein H. Membrane Innovation in Dialysis. *Contrib Nephrol*. 2017;191:100-114.
44. Lorenzin A, Neri M, Lupi A, et al. Quantification of internal filtration in hollow fiber hemodialyzers with medium cut-off membrane. *Blood Purification*. 2018;46(3):196-204.
45. Zuo L, Wu J, Yao Q, et al. Classification of hemodialyzer clinical performance. Zhongguancun Nephropathy Blood Purification Innovation Alliance. 2021:1-4.
46. Clinical Data. Baxter Renal Care. <https://renalcare.baxter.com/data-insights/clinical-data>. Accessed April 26, 2022.
47. Kandi M, Brignardello-Petersen R, Couban R, Wu C, Nesrallah G. Clinical Outcomes with Medium Cut-Off Versus High-Flux Hemodialysis Membranes: A Systematic Review and Meta-Analysis. *Can J Kidney Health Dis*. 2022;9:1-16. DOI: 10.1177/20543581211067087
48. Ronco C, La Manna G. Expanded Hemodialysis: A New Therapy for a New Class of Membranes. *Contrib Nephrol*. 2017;190:124-133. doi:10.1159/000468959
49. Wolley M, Jardine M, Hutchison CA. Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol*. 2018;13(5):805-814. doi:10.2215/CJN.10110917
50. EUTOX Uremic Solutes Database. June 2022. [Uremic-toxins.org](https://www.uremic-toxins.org)

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Baxter International Inc. oder ihrer Tochtergesellschaften.
Gambro Dialysatoren GmbH
**Baxter Deutschland GmbH**